

Histology Outlet Mall Enterprises, Inc.

P.O. Box 407

Auburn, IN 46706

800-474-9477

FAX 260-333-3017

IMMUNOPEROXIDASE CONTROLS

Type	CAT #	Box of 10	Box of 20	Box of 100
		A	B	C
ALPHA-1-FETOPROTEIN (AFP)	1055	\$41	\$62	\$243
B-CELL, CD45RA	1080	\$41	\$62	\$243
BCL-2	1085	\$41	\$62	\$243
BCL-6	1095	\$41	\$62	\$243
BROMODEOXYURIDINE (BrdUrd)	1090	\$41	\$62	\$243
CD 3	1103	\$41	\$62	\$243
CD 5	1105	\$41	\$62	\$243
CD 8	1108	\$41	\$62	\$243
CD 10	1110	\$41	\$62	\$243
CD 20 (L26)	1120	\$41	\$62	\$243
CD 30 (Ki-1)	1130	\$41	\$62	\$243
CD 43	1143	\$41	\$62	\$243
CD 45RA, B-CELL	1080	\$41	\$62	\$243
CDX2	1102	\$41	\$62	\$243
CALRETININ	1150	\$41	\$62	\$243
CARCINOEMBRYONIC ANTIGEN (CEA)	1160	\$41	\$62	\$243
CHROMOGRANIN	1190	\$41	\$62	\$243
CYTOKERATIN-MONOCLONAL	1235	\$41	\$62	\$243
CYTOKERATIN-POLYCLONAL	1236	\$41	\$62	\$243
CYTOMEGALOVIRUS (CMV)	1240	\$41	\$62	\$243
DESMIN	1255	\$41	\$62	\$243
EPITHELIAL MEMBRANE ANTIGEN (EMA)	1285	\$41	\$62	\$243
FACTOR 8	1290	\$41	\$62	\$243
FACTOR 13	1295	\$41	\$62	\$243
GLUCAGON	1310	\$41	\$62	\$243
HEPATITIS B CORE ANTIGEN	1360	\$41	\$62	\$243
HEPATITIS B SURFACE ANTIGEN	1365	\$41	\$62	\$243
HMB 45	1370	\$41	\$62	N/A
HERPES SIMPLEX VIRUS, TYPE I (HSV-1)	1380	\$41	\$62	\$243
HERPES SIMPLEX VIRUS, TYPE II (HSV-2)	1385	\$41	\$62	\$243
INSULIN	1420	\$41	\$62	\$243
KAPPA	1500	\$41	\$62	\$243
LAMBDA	1540	\$41	\$62	\$243
LEUCOCYTE COMMON ANTIGEN, CD45 (LCA)	1545	\$41	\$62	\$243
LN 2	1560	\$41	\$62	\$243
NEURON SPECIFIC ENOLASE (NSE)	1610	\$41	\$62	\$243
p53 PROTEIN	1640	\$41	\$62	\$243
PROSTATE SPECIFIC ANTIGEN (PSA)	1680	\$41	\$62	\$243
PROSTATIC ACID PHOSPHATASE (PSAP)	1685	\$41	\$62	\$243
S-100	1740	\$41	\$62	\$243
SOMATOSTATIN	1760	\$41	\$62	\$243
SYNAPTOPHYSIN	1770	\$41	\$62	\$243
T-CELL, CD45RO	1780	\$41	\$62	\$243
VIMENTIN	1870	\$41	\$62	\$243

Every box of 10 contains 10 unstained serial sections and one stained slide.

Every box of 20 contains 20 unstained serial sections and one stained slide.

Every box of 100 contains 100 unstained serial sections and one stained slide.

In addition, each box comes with the staining procedure inside.

All Histopathology controls are on aminosilane (positive charged) slides.

***4% additional charge will be applied if paying by credit card**

***If you use a Ventana staining machine, we recommend requesting the M1000 slides - there is a \$10 upcharge.
Place "M" in front of catalog number*

***If you prefer tissue at the top of the slides, no problem!! Lead time may be extended, but no extra charge! Follow the catalog number by "-T"*